



ABI Genotyping total solution


SNP Genotyping Technologies Available from AB

SNP ANALYSIS

PLATFORM


Chromosome & Genome-Wide LD Studies

- SNIPlex™ on electrophoretic genetic analyzer
- High throughput




Candidate Gene / Regional LD Studies

- TaqMan® on Sequence Detection System
- Cost-effective, easy to use, validated assays



Dense Mutational Spectra

- Resequencing Primers on electrophoretic genetic analyzer
- Dual use: mutation discovery & assay






ABI Genotyping platform

- 5' Nuclease Assay = TaqMan Probe mismatch
 - Allele Discrimination
 - 7900HT, 7700, 7000
- Sequencing = actually determine/compare sequences
 - 3730, 3100, 310
 - VariantSeqr resequencing primer
 - Sequencing reagent:
 - SeqScape v2.1
- Multiplexed OLA/PCR Assay Extension = SNPlex assay
 - 3730, (3100)
 - GeneMapper v 3.5



ABI Tool Box of Technologies...

Sequence Detection Systems



5' Nuclease Assay

Electrophoresis Platforms




***Resequencing
SNPlex Assay***

TaqMan 5' Nuclease Assay

For Allele 1 - only VIC™ dye signal is generated




For Allele 2 - only FAM™ dye signal is generated



Detection of fluorescent dye signal
by sequence detection system



Allelic Discrimination Cluster Data View





TaqMan[®] Assay → Simple Workflow

One step, single tube assay

Setup

Pipette: Master Mix + Assay mix + DNA Sample & Seal plate



Cycle




Read



**ABI PRISM[®] 7900HT SDS
with integrated robotics**



Analyze



Three Solutions for SNP design



- Ready made assays
- 網路選購
- Updated Human Genome information



- 代客設計與測試 for SNPs
- 85% 成功率



- Design primer or probe by software



Applied Biosystems


針對人類3萬個基因,推出20萬個SNP-AoD基因分析試劑組

SNP Assays

Assays-on-Demand products

Gene Expression


Assays-on-Demand products



SNP
Gene Expression
Assays-by-DesignSM Service



AB SNP Validation Overview





SNP AoD Selection Goals

- A genome-wide SNP marker set useful for LD mapping
 - 200,000 5' Nuclease SNP Assays
 - SNPs in gene-regions / gene centric approach
 - Minor allele frequency >10% in one population
 - Coverage of 1 SNP per 10 Kb of gene length



AbD Position as


- Custom Assay
 - Unlike off the shelf you get what you want
- Convenience
 - Unlike oligos in many tubes
 - All in one tube
 - Assay is designed and optimized (standard thermal cycling conditions)
- High Quality Control
 - Mass spec
 - HPLC purified probe
 - **Functional tested??** •85% 成功率



Do-It-Yourself vs. The Assays-by-DesignSM Service

Streamlining Your Workflow

Do it yourself



Assays-by-Design Service



* Not required when using standard SDS protocols





Using SNPbrowser™ to Help Design Genetic Association Studies


Selecting SNPs for Association Studies

Factors to be considered when selecting SNPs:

- minor allele frequency
- distance between adjacent SNPs
- location with respect to exons or regulatory elements
- non-coding or coding
- local LD architecture

View Chr. Position Bases: 72,143,471 Distance LDUs: 1260.399

HUGO Name Search * Batch Search Order SNPs Print



View Chr.

Position

Distance

12

Bases: 72,052,308

LDUs: 1200.271

2.051



Zoom Out

HUGO Name

tph2*



Search *



Batch Search



Order SNPs



Print

Cauc. LD B

SNP Details

SNP Name

hCV245410

Location

Chr. 12 Base Pos. 72,053,260

Min. Allele Frequencies

AfAm: 0.20

Cauc: 0.43

Chin: 0.41

Japn: 0.48

Chromosome 12

Ave. Power Scale



hCV22272411

hCV2904133

hCV12057924

hCV245410

hCV8872342

hCV15836061

hCV8376173

hCV8376146

hCV8872308

hCV8376042

hCV8872233

hCV11407441

(58021111)98249000

We can deselect redundant SNPs based on:

- LDU map information
- Minor allele frequency data

72,124,598

72,166,185

Bases

1199.501

1201.687

Caucasian LDUs



Using SNPbrowser to Help Design Genetic Association Studies

Which SNPs Should be Chosen for Genotyping?

For reasons of efficiency, we may not want to type all 10 of these SNPs (at least in the first round of genotyping).

Using allele frequency information and the LDU scale we can deselect those SNPs that might not give much additional information.

	hCV12057924	0.13		
	hCV245410	0.43		
	hCV8872342	0.08		
	hCV15836061	0.34		
	hCV8376173	0.35		
	hCV8376146	0.33	}	these SNPs have similar minor allele frequencies and are coincident on the LDU map - probably don't all need to be genotyped
Deselected	hCV8872308	0.08		
assays	hCV8376042	0.22	}	
	hCV8872233	0.24		
	hCV11407441	0.23		

similar minor allele frequency to adjacent SNP, but separated by small LDU distance



ABI Genotyping platform


- 5' Nuclease Assay = TaqMan Probe mismatch
 - Allele Discrimination
 - 7900HT, 7700, 7000
- Sequencing = actually determine/compare sequences
 - 3730, 3100, 310
 - VariantSeqr resequencing primer
 - Sequencing reagent:
 - SeqScape v2.1
- Multiplexed OLA/PCR Assay Extension = SNPlex assay
 - 3730, (3100)
 - GeneMapper v 3.5



When is sequencing the best technology for detecting variants?

- Discovering new mutations
- Detection of insertions/deletions as well as novel SNPs
- Dense mutational spectrum
- Evaluating known variants

VariantSEQr™ Resequencing System Complete Gene Primer Sets



- Includes all Genbank RefSeqs (transcripts) for each gene → ‘supertranscript’
- Primers for all exons (coding and non-coding)
- ‘Regulatory region’ = 1kb 5’ of exon 1 start
- Splice junctions
- Primers tailed with M13 universal F and R sequencing primers
- Annotation of regions where coverage is < 100%



An Application Product That Is Applied Biosystems Optimized



AmpliTaq Gold[®]
DNA Polymerase



✓ *Resequencing Primers:*

- ❖ *Ready-To-Pipet Into PCR Reaction*
- ❖ *Enough Liquid Volume For 500 DNA Samples*

✓ *Protocol:*

- ❖ *PCR Thermal Cycling Profile*
- ❖ *PCR Cleanup*
- ❖ *Sequencing Thermal Cycling Profile*
- ❖ *Sequencing Cleanup*

✓ *Optimized For Use On:*

- ❖ *Applied Biosystems 9700 thermal cycler*
- ❖ *AmpliTaq Gold[®] DNA Polymerase*
- ❖ *BigDye[®] Terminator v3.1*
- ❖ *Applied Biosystems 3730 & 3730xl & 3100 systems*
- ❖ *Applied Biosystems SeqScape[®] Software 2.1*



VariantSEQr™ System Summary

- Validated PCR primers to provide consistent and specific amplification
- Primers pre-mixed, diluted and aliquoted for easy use
- Streamlined and optimized protocols for efficient workflow
- Automated and highly accurate data analysis



ABI Genotyping platform

- 5' Nuclease Assay = TaqMan Probe mismatch
 - Allele Discrimination
 - 7900HT, 7700, 7000
- Sequencing = actually determine/compare sequences
 - 3730, 3100, 310
 - VariantSEQr resequencing primer
 - Sequencing reagent:
 - SeqScape v2.1
- Multiplexed OLA/PCR Assay Extension = SNPlex assay
 - 3730, (3100)
 - GeneMapper v 3.5



Complete Solution for High Throughput SNP Genotyping

Locus Specific Assays

Fixed Content

Human Linkage Mapping Set


Custom Content

Instrumentation

Applied Biosystems 3730/3730xl

POP-7™, 36 cm array

15 minutes per run



Universal Reagents

OLA/PCR modules

Reporter Probes

Multi-color fluorescence

Mobility Modifiers

Data Analysis

Data Collection Software v2.0

GeneMapper™ Software v 3.5



Assay Products

- **SNIPlex™ System Custom: (Flexible, made-to-order):**
 - Customer specified/selected
 - AB designs probes & multiplex groups, synthesizes oligos, pools probes, ships to customer
- **SNIPlex™ System Fixed (off-the-shelf):**
 - Human Linkage Mapping Set ~ 3500 markers
 - Additional Human and Mouse Sets planned



Basic Assay Steps Involved in SNPLex™ System

Ligation

Generation of genotype (GT) specific products through multiplex oligonucleotide ligation reaction (OLA)



Amplification


Multiplex PCR with universal primers



Detection

Hybridization of universal ZipChute™ probes to amplicons and identification of eluted ZipChutes by CE

Universal ZipChute™ Probes



- The same **Universal ZipChute** mix is used for hybridization to each multiplex reaction, it is not SNP specific
- ZipChute probes that specifically hybridize to genotype specific amplicons are eluted and analyzed by capillary electrophoresis
- ZipChute identification is used for genotype determination



Whole Product View

- **OLA/PCR Reagent Kits (Universal)**
 1. Phosphorylation Module
 2. Ligation Module
 3. Amplification Module
 4. Purification Module
 5. Hybridization Module
 6. Controls
 7. Standards Module
- **SNP specific probes (Assays)**
- **GeneMapper™ 3.5, Data Collection 2.0, BioTrekker™ 1.0 (optional)**
- **Electrophoresis instruments & consumables**

SNIPlex™ System OLA/PCR Assay

Universal PCR Priming site

GER

ZipCode₁

A

ASO_x

G

ASO_y

ZipCode₂

gDNA
Target

C

GER = Genome Equivalent Region
ASO = Allele Specific Oligo
LSO = Locus Specific Oligo

Universal PCR Priming site

P

LSO

GER

1. OLA

P

A

G

Ligation Product Formed

C

(Homozygote shown in this case)

2. Clean-up



PCR & ZipChute™ Probe Hybridization

3. Multiplexed Universal PCR

Univ. PCR Primer




4. Capture double stranded DNA



5. Denature double stranded DNA


6. Wash away one strand

7. ZipChute™ Probe Hybridization



Detection

8. Wash and release ZipChute™ Probes




9. Load on CE instrument




2004: Experimental Design



**Applied Biosystems
3730xl DNA Analyzer**



**SNIPlex™ Genotyping
System**



CELERA DISCOVERY SYSTEM™
The essential tool for the life science researcher




**ABI PRISM® 7900HT
Sequence Detection System**



**TaqMan® SNP Genotyping Assays
(Formerly Assays-on-Demand™
products)**



Disease Gene Hunting:



AB offer total solution



以計畫的需求來決定適用方案： 唯有ABI能提供**SNP Total Solutions**

	# SNP Markers	# Samples	Recommended Platform
Whole Genome	1000's-500K	1000's	SNPlex
Linkage Mapping	1000's	100's-1000	SNPlex
Candidate Gene/Region	10-1000's	500-1000's	TaqMan/SNPlex
ADME/Tox	100's	10-1000's	TaqMan



SNPlex vs TaqMan

	SNPlex	TaqMan
實驗流程	Ligation Probe Hybridization Detection	Hybridization Detection
Throughput (genotypes per day)	437,760 - 1,300,000	27,648 - 255,000
Running cost (per SNP)	0.02-0.08	0.4-0.5



Why customers will not switch from Taqman to SNPLex

- Don't have CE instr. or budget
- Don't have robotics or budget
- Have projects with critical SNPs (i.e. can't substitute)
- Have projects with small numbers of SNPs (< 50)
- Working exclusively with non-Human species
 - SNPLex can support this, but optimized for Human



Why Customers Might Switch from Taqman to SNPlex...

- Attractiveness of price/GT
 - This can be the wrong reason based upon approach/study needs
- Already has CE platform
- Wants to do sequencing, re-sequencing, linkage or other fragment analysis applications
- Existing expertise/hardware for automation
- Project size is large
 - Dictates need for high throughput capability

Two Technology Choices

Choose TaqMan


1. To take advantage of **singleplex** technology for “mix & match” of individual SNPs
2. When you are **limited to a small** number of SNPs to genotype
3. For simplified **one-step workflow** and **convenience** of validated assays
4. Additionally require **gene expression** capability in your laboratory

Choose SNPlex

1. To take advantage of **multiplex** technology to lower overall project cost
2. When you have desired “set” of 48 or more SNPs to genotype
3. Use **robotics/automation** in your laboratory
4. Additionally require **sequencing** capability in your laboratory

Complete SNP genotyping product portfolio:
Candidate gene/candidate region studies
Whole chromosome/genome studies

SNIPlex™ Assay Workflow Overview



5' Nuclease based solution is simplest to implement

- Easiest to implement
 - Assays by DesignSM service
 - Ready to use Assays on DemandTM
 - Universal master mix
 - Minimum of pipetting steps
 - no transfer or additions post PCR.
- Easiest to transfer between collaborators
 - PCR infrastructure + 7900, 7700, or 7000
 - No multiplexing required

Assays are portable
from study to study,
lab to lab



AB SNPS- 100% Ready To Use

- 200,000
 - fully-validated
 - true SNPs
 - Real allele frequency data
 - Hardy-Weinberg Equilibrium
 - Gene coverage plus
 - Ready-to-use assays

Choose it and use it!




The Strategy


- Shift customer focus from per genotype costs to project costs.
- AB providing enabling resources that are redefining project time and cost




ABI 基因體研究相關的全貌




GENE EXPRESSION



SNP GENOTYPING



SEQUENCING



RESEQUENCING



CELERA DISCOVERY SYSTEM
The essential tool for the life science researcher

Are You a Current CDS Subscriber? What's New Publications Service Offerings

Market and Distributed by **ABI Applied Biosystems**